


PLEASE READ THESE INSTALLATION INSTRUCTIONS THOROUGHLY PRIOR TO INSTALLING THE SEAT

THIS PRODUCT SHOULD BE INSTALLED BY QUALIFIED PERSONNEL

1. Determine desired TUB BENCH location.

The Bench must be located on the bathtub so that the swing-down legs are resting on the flat of the tub floor, and that the adjustable leveling feet are positioned on a clean, flat surface.



Remove the Tub Bench from the shipping carton.

Lay the Bench upside down on a clean flat surface.

Rotate the Swing Down Legs so they are positioned perpendicular to the seat top.

Line up the hole in the Swing Down Leg tubing with the hole in the center of the Swing Down Leg Bracket.

Place the threaded Thumb Screw into the center hole on the Swing Down Leg Bracket and turn clockwise until "hand tight" to secure the Swing Down Legs in a position perpendicular to the seat top.

2. Determine Bench height from the floor and location relative to the tub and bathroom finished floor.


ADA regulations specify a finished height of 17" to 19" above the bathroom finished floor.

Holding the Bench in its horizontal position in the tub, adjust the Tub Clamps so they will fit over the front of the tub ledge.

The length of the "Swing Down Legs" may require adjusting so the Bench rests in a flat horizontal position. To change the height of the legs, adjust the length of the adjustable leveling feet and tighten the lock nut to secure them in place.

Tighten the Tub Clamp to secure the Bench on the tube ledge.

Do not "over-tighten" the Tub Clamp, which may distort the Tub Clamp Bracket.



Make certain that the Leveling Feet are resting on a clean, FLAT SURFACE ON THE TUB FLOOR.

3. Check Seat for:

- 1) Correct finished height.
- 2) Tub Clamp properly secured on tub ledge.
- 3) "Swing-Down Legs" locked in down position @ 90 degrees, perpendicular to the tub floor.
- 4) The adjustable leveling feet are positioned on a clean, flat surface on the tub floor, not resting on the curved/radius corner of the tub floor.

Once secured to the tub, make certain the seat does not shift from its set position in any direction, before use.

THIS COMPLETES THE INSTALLATION OF THE TUB SEAT.

Removing the Bench:

When removing the Bench, reverse the above procedures:

1. Loosen (turn counter-clockwise) the Tub Clamp to remove the seat from the tub ledge.
2. Lay the Bench on a clean, flat surface so the Swing Down Leg Bracket is accessible.
3. Loosen (turn counter-clockwise) the threaded Thumb Screw to remove it from the center hole of the Swing Down Leg Bracket.
4. Rotate the Swing Down Legs until they rest on the underside of the Bench top.
5. Store this product in a clean, dry location for future use.

For technical assistance, contact WingIts....