



PROFECT® 64 HDQC2®

Daily Cleaner Disinfectant



Profect® 64 hdqc2 is a one step¹ hospital disinfectant cleaner¹ formulated to kill a broad spectrum of microorganisms such as HIV-1 (AIDS Virus), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV) on hard, nonporous surfaces. With a neutral pH, Profect 64 hdqc2 disinfects without damaging or dulling floor finish even with repeated use. Available in a refreshing citrus scent. Profect 64 hdqc2 is cost-effective and reduces restocking intervals for routine high touch surface cleaning. Use Profect 64 hdqc2 daily as part of a simple and effective cleaning and disinfection¹ program for your entire facility.

EPA Registration Number: 5741-35

PROFECT® 64 HDQC2®

FEATURES

- pH neutral formulation designed to provide effective cleaning, deodorizing, and disinfection for all hard, nonporous environmental surfaces
- Bactericide/Virucide*/Mildewstat
- 2 minute sanitizing³ activity***
- Sanitizes soft surfaces^{2,3} in 10 minutes

BENEFITS

- Economical to use - dilute at 2 ounces per gallon of water
- Meets surface disinfection recommendations of OSHA's Bloodborne Pathogens Standard
- One product for all your cleaning and disinfection needs
- Kills SARS-CoV2⁸ the virus that causes COVID-19

¹When used according to disinfection directions for all listed organisms

²100% cotton and 100% polyester

³When used according to the soft surface sanitizer directions for spot treatment of precleaned soft surfaces

⁸On hard, non-porous surfaces

^{*}When used as directed for sanitizing hard, non-porous, non-food contact surfaces

^{*}Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), HIV-1 (AIDS Virus), Influenza A Virus, and Influenza B Virus

^{***}*Staphylococcus aureus* (ATCC 13048) and *Enterobacter aerogenes* (ATCC 6538)

EFFICACY DATA

FOR USE AS A ONE-STEP¹ DISINFECTANT, DEODORIZER, AND CLEANER ON HARD, NON-POROUS SURFACES

When used as directed, this product is effective against *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Bordetella bronchiseptica*, *Enterobacter aerogenes*, *Escherichia coli* O157:H7, Extended-Spectrum beta-lactamase (ESBL) producing *Escherichia coli*, *Klebsiella pneumoniae* (ESBL), *Klebsiella pneumoniae* (NDM-1), *Listeria monocytogenes*, Methicillin Resistant *Staphylococcus aureus* (MRSA), Multi-drug Resistant (MDR) *Acinetobacter baumannii*, *Neisseria gonorrhoeae*, *Proteus mirabilis*, *Proteus vulgaris*, *Salmonella enterica*, *Salmonella enterica* subspecies enterica serovar Pullorum, *Salmonella enterica* subspecies enterica serovar Typhimurium, *Serratia marcescens*, *Shigella flexneri* serotype 2B, *Shigella sonnei*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, Vancomycin Resistant *Enterococcus faecalis* (VRE), *Vibrio cholerae*, and *Yersinia enterocolitica*.

*VIRUCIDAL ACTIVITY ON HARD, NON-POROUS SURFACES

When used as directed, this product is effective against Hepatitis C Virus (HCV), Hepatitis B Virus (HBV), HIV-1 (AIDS Virus), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Coronavirus, Human Rotavirus, Influenza A Virus, Influenza B Virus, SARS-CoV-2 and Vaccinia virus.

TO CONTROL MILDEW

When used as directed this product is an effective mildewstat against *Aspergillus niger*.

SPECIFICATION DATA

Dilution	1:64
pH	6.0 – 9.0
Color	Red
Fragrance	Refreshing Herbal
Specific Gravity @ 24°C/75°F	0.999
Flash Point	None
Stability: Shelf @ 24°C/75°F	1-year minimum
Stability: Freeze/Thaw	3 cycles
EPA Registration Number	5741-35

ORDERING INFORMATION



ADDITIONAL RESOURCES

Scan the QR code to access this product's safety information and documentation.



SAFETY: For institutional and industrial use only. Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at www.spartanchemical.com. Workplace labels are available in English, Spanish, and French.

GUARANTEE: Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.