

KTC1 WET END PARTS


The following part numbers are the pump models with liquid end code KTC1 in the model  
EXAMPLE: LB64SA-KTC1-XXX

Use the fourth digit in the model number from your pump to identify the part numbers below

EXAMPLE MODEL:

|          |   |   |   |   |   |   |   |   |   |   |    |   |    |    |    |
|----------|---|---|---|---|---|---|---|---|---|---|----|---|----|----|----|
| Position | 1 | 2 | 3 | 4 | 5 | 6 | - | 7 | 8 | 9 | 10 | - | 11 | 12 | 13 |
| Model:   | L | B | 6 | 4 | S | A | - | K | T | C | 1  | - | X  | X  | X  |

For example the pump head (ITEM #1) for this pump would be:  
Part Number: L0200400-PVD



| PART NUMBER                                                                      | QTY | DIGIT IN POSITION NUMBER<br>FOUR OF THE MODEL NUMBER       |
|----------------------------------------------------------------------------------|-----|------------------------------------------------------------|
| ITEM NO. 1: PUMP HEAD (PVD)                                                      |     |                                                            |
| L0200900-PVD                                                                     | 1   | 2                                                          |
| L0200300-PVD                                                                     | 1   | 3                                                          |
| L0200400-PVD                                                                     | 1   | 4                                                          |
| ITEM NO. 2: DIAPHRAGM                                                            |     |                                                            |
| L0300900-THY                                                                     | 1   | 2                                                          |
| L0301000-THY                                                                     | 1   | 3                                                          |
| L0301100-THY                                                                     | 1   | 4                                                          |
| ITEM NO. 3: DEFLECTION PLATE                                                     |     |                                                            |
| L2100200-FPP                                                                     | 1   | 2                                                          |
| L2100300-FPP                                                                     | 1   | 3                                                          |
| L2100400-FPP                                                                     | 1   | 4                                                          |
| ITEM NO. 5: SHIM (0-2 REQUIRED)                                                  |     |                                                            |
| L9901200-BRS                                                                     | 0-2 | ANY DIGIT                                                  |
| ITEM NO. 7 & 8: PUMP HEAD SCREWS AND WASHERS<br>KIT (CONTAINS BOLTS AND WASHERS) |     |                                                            |
| L3900200-000                                                                     | 1   | 2, 3 OR 4                                                  |
| ITEM NO. 9: SUCTION VALVE                                                        |     |                                                            |
| L3101TC1-PVD                                                                     | 1   | ANY DIGIT                                                  |
| ITEM NO. 10: DISCHARGE VALVE                                                     |     |                                                            |
| L3201TC1-PVD                                                                     | 1   | ANY DIGIT                                                  |
| ITEM NO. 11: BLEED VALVE                                                         |     |                                                            |
| L3300T01-PVD                                                                     | 1   | ANY DIGIT                                                  |
| ITEM NO. 12: FOOT VALVE/STRAINER                                                 |     |                                                            |
| J60716                                                                           | 1   | ANY DIGIT                                                  |
| ITEM NO. 13: INJECTION VALVE                                                     |     |                                                            |
| J61020                                                                           | 1   | ANY DIGIT                                                  |
| REPLACEMENT FUSES                                                                |     |                                                            |
| L9707300-000                                                                     | 1   | E+ SERIES AND E SERIES PUMPS<br>WITH FUSE HOLDER           |
| L9711200-000                                                                     | 1   | A+, C, C+ AND E SERIES PUMPS<br>WITH PC BOARD MOUNTED FUSE |
| L9704700-000                                                                     | 1   | SERIES E-DC PUMPS                                          |
| L9700700-250                                                                     | 1   | SERIES MP PUMP-CIRCUIT                                     |

There is a KOPkit (rebuild kit) available for your pump that includes items  
1, 2, 5, 7, 8, 9, 10 and a bleed valve o-ring

The part number for the KOPkit can be found on your pump data label or  
you can determine the KOPkit part number from your pump model number  
by using the fourth, seventh, eighth, ninth and tenth positions for the model string

EXAMPLE:

|               |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |          |
|---------------|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|----------|
| Model         | L | B | 6 | 4 | S | A | - | K | T | C | 1 | - | X | X | X |          |
| KOPkit Number |   |   |   | 4 |   |   |   | K | T | C | 1 |   |   |   |   | = K4KTC1 |