



Think ahead.

# Tork Advanced Facial Tissue Flat Box

TF6810



## Description

Tork Advanced Facial Tissue is an extra soft and absorbent 2-ply facial tissue which offers the right combination of performance and value. The flat box style is ideal for the vanity or in-wall dispensers.

- Soft protective poly shield opening protects the tissues and provides easy, one-at-a-time dispensing while reducing lint.
- Environmental Label
- Advanced
- One hygienic sheet at a time
- 100 % Recycled
- Natural Color Fibers & Dye Free

## Product Certifications

SmartOne®



Tork Advanced

## Product Details

System	F1
--------	----

## Shipping Data

	Consumer Units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	7016646002486	10073286614886	7322540850581
Packaging Material	Box	Carton	-
Pieces	100	30 (30 CON)	243000 (81 TRP)
Height	2.09 in	10.75 in	96.75 in
Length	4.53 in	18 in	54 in
Width	8.46 in	14.25 in	42.75 in
Gross Weight	0.37 lb	11.47 lb	928.59 lb
Net Weight	0.26 lb	7.89 lb	639.33 lb
Volume	0.05 ft	1.59 ft	129.25 ft
Layers Per Pallet	-	-	9
TRP Per Layer	-	-	9

## Environmental Information

### Content

The product is made from

Recycled fibers  
Chemicals

The packaging material is made from paper or plastic.

### Material

Recycled fibers

Recycling of paper is an efficient use of resources as the wood fibers are used more than once.

High demands are put on quality and purity of recovered fibers, considering each step of the chain (collecting, sorting, transporting, storage, use), to ensure safe and hygienic products.

Recovered paper can be produced both from collected newsprint, magazines and office waste. The choice of recovered paper grades, is made for each product, depending on its specific requirements on performance properties and brightness. The paper is dissolved in water, washed and treated with chemicals under high temperature and screened to separate out impurities.

Bleaching is a cleaning process of the fibers that is often used. The aim is then to achieve a bright pulp, but also to get a certain purity of the fiber in order to achieve the demands for hygiene products and in some cases to meet the requirements for food safety.

Bleaching of the recovered pulp is made with chlorine-free bleaching agents (hydrogene peroxide and sodium dithionite). Except for Natural Napkins that are unbleached.

For bleached products we use bleaching agents (to increase the brightness of pulp from recovered paper).

### Chemicals

All chemicals (process aids as well as additives) are assessed from an environmental, occupational health and safety and product safety point of view.

To control product performance we use additives:

- Wet strength agents (for Wipers and Hand Towels)
- Dry strength agents (is used together with mechanical treatment of the pulp to make strong products like wipers)
- For colored papers dyes and fixatives (to secure perfect fastness of the color) are added
- For printed products printing inks (pigments with carriers and fixatives) are applied
- For multi ply products we often use water soluble glue to secure the integrity of the product

In most of our mills we do not add optical brighteners but it often occurs in recovered paper since it is used in printing paper.

We do not use softeners for professional hygiene products.

High product quality is secured through quality and hygiene management systems throughout production, storage and transport.

In order to maintain a stable process and product quality the paper manufacturing process is supported by the following chemicals/ process aids:

- defoamers (surfactants and dispersing agents)
- pH-control (sodium hydroxide and sulphuric acid)
- retention aids (chemicals that help to agglomerate small fibers to prevent fiber loss)

# Tork Advanced Facial Tissue Flat Box

TF6810

- Coating chemicals (that help to control the creping of the paper to make it soft and absorbent)

To reuse broke and to utilize recovered fibers we use:

- Pulping aid (chemicals that help to repulp wet strong paper)
- Flocculation chemicals (that help to clean out printing inks and fillers from recovered paper)
- Bleaching agents (to increase the brightness of pulp from recovered paper)

In the cleaning of our waste water we use flocculation agents and nutrients for the biological treatment to secure that no negative impact on water quality comes from our mills.

<b>Environmental certification</b>	This product is certified for FSC®.
<b>Packaging</b>	Fulfilment of Packaging and Packaging Waste Directive (94/62/EC): Yes
<b>Article creation date and latest article revision</b>	Date of issue: 19-04-2019 Revision date: 02-04-2024
<b>Production</b>	This product is produced at Essity Professional Hygiene NA Blending - US mill.
<b>Destruction</b>	This product is mainly used for personal hygiene and can be collected together with household waste.

**Essity North America Inc., Cira Centre, Suite 2600 2929 Arch Street, Philadelphia, PA 19104, USA**