

Aquanomic BioCare Solid Detergent

CLASSIFICATION

Concentrated laundry detergent

PRODUCT DESCRIPTION

This enzymatic laundry detergent is designed to combat tough food and body soils, resulting in clean, odor-free linens. Washes at a lower pH than traditional alkaline detergents, reducing carryover and allowing premium bleach performance. Can be used in a combined suds/bleach step for increased enzyme contact time.

USE

FOR INDUSTRIAL USE ONLY – NOT FOR HOUSEHOLD USE

DIRECTIONS FOR USE

1. Remove product container cap. Use container handle. Place container open end down in the Ecolab Solid Laundry Dispenser.
2. To remove used container, carefully lift from dispenser. Turn container upright.
3. Use cap from new container to close the used container. Do not reuse empty container.
4. If product contacts hands, wash off immediately. Use product for intended purpose or dispose in accordance with regulations. Rinse empty container with water 3 times and recycle empty container.

Temperature: Recommended wash temperature 100°F - 160°F

Use only in an Ecolab dispenser.

Consult your Ecolab Specialist or call 1-800-35 CLEAN (352-5326) for specific use directions, service, or additional information.

PHYSICAL & CHEMICAL PROPERTIES

Appearance: Cast Solid
Color: Opaque, white
Odor: Not significant
pH: 10.0 – 12.0, 1%

Ingredients:

Sodium carbonate	497-19-8
C12-C14 pareth-7	68551-12-2
Sodium chloride	7647-14-5
Sodium bicarbonate	144-55-8
Poly(oxy-1,2-ethanediyl), .alpha.-isotridecyl-.omega.-hydroxy-	9043-30-5
Poly(oxy-1,2-ethanediyl), a-(2-propylheptyl)-w-hydroxy-	160875-66-1
Polyethylene glycol	25322-68-3
Ethylenediamine tetraacetate	64-02-8
Dodecylbenzenesulfonic acid, sodium salt	25155-30-0

%P: Formula ingredients contain no phosphorus

PACKAGING & PRODUCT NUMBERS

Size	Product Number
4 - 3.6 KG (8 LB)	6101752